

Liste des produit QuikRead go Instrument

• Panels de tests :

- o **iFOBT** : Détection et quantification de sang occulte dans les selles, Le système affiche des valeurs quantitatives de 75 à 1000 ng/ml dans la solution d'analyse. Les valeurs quantitatives peuvent aussi être affichées sous la forme suivant $\mu\text{g/g}$ de selles.
- o **StrepA** : Détection de Streptococcus pyogenes dans le frottis de gorge
- o **wR CRP** : WIDE Rang CRP, L'intervalle de mesure du test de la CRP s'étend de 0.5 à 300 mg/l
- o **wCRP+Hb** : L'intervalle de mesure de la CRP s'étend de 5 à 200 mg/l et celui de l'hémoglobine s'étend quant à lui de 50 à 245 g/l (3,1 à 15,2 mmol/l).

QuikRead go wrCRP	Cat.No.
QuikRead go wrCRP, 50 tests including capillaries (10 μl) and plungers	146521
QuikRead go wrCRP, 500 tests	146522
QuikRead CRP Control	68296
QuikRead go CRP Control High	137071
Ultrasensitive CRP Control	68257
QuikRead go wrCRP+Hb	Cat.No.
QuikRead go wrCRP+Hb, 50 tests including capillaries (10 μl) and plungers	146255
QuikRead CRP Control	68296
QuikRead go CRP Control High	137071
Ultrasensitive CRP Control	68257
QuikRead go Hb Control	141154
QuikRead go Strep A	Cat.No.
QuikRead go Strep A, 50 tests including controls	135883
QuikRead go iFOBT	Cat.No.
QuikRead go iFOBT, 50 tests	136189
QuikRead FOB Sampling Set, 50 tests	138822
QuikRead FOB Control Quantitative	144675
QuikRead FOB Positive Control	06027

